

Physiological, Microscopic, and Molecular Approaches to Characterize Methane-Oxidizing *Methylosinus* sp. Strain C49: Culture Revival, Fluorescence Imaging Optimization, and Purity Assessment

Jameil Magomnang Jr.¹, Akihiko Terada^{1,2,*}

¹Department of Applied Physics and Chemical Engineering, Graduate School of Engineering, Tokyo University of Agriculture and Technology, 2-24-16 Naka, Koganei, 184-8588, Tokyo, Japan

²Institute of Global Innovation Research, Tokyo University of Agriculture and Technology, 3-8-1 Harumi, Fuchu, 185-8538, Tokyo, Japan

*Author to whom correspondence should be addressed:
E-mail: akte@cc.tuat.ac.jp

(Received April 06, 2026; Revised June 30, 2026; Accepted September 16, 2206)

Supplementary Material




Fig. S1: Schematic of the continuous optical density (OD₆₀₀) monitoring system showing the culture tube, non-contact turbidimeter, OD-control unit, and PC running OD Viewer software




Fig. S2: Shaker-mounted cultivation system equipped with the non-contact turbidimeter for continuous OD₆₀₀ monitoring of MC49 cultures




Fig. S3: Visual appearance of liquid cultures of *Methylosinus* sp. strain C49 in nitrate mineral salts (NMS) medium under methane-rich conditions. Panel (A) shows the negative control vial without inoculation, displaying a clear, transparent medium indicative of no microbial growth. Panel (B) depicts the culture immediately after inoculation with an initial OD₆₀₀ of approximately 0.32, reflecting the initial cell density. Panel (C) shows the culture after 14 days of incubation at 30°C with shaking at 130 rpm, exhibiting noticeably increased turbidity corresponding to a maximum OD₆₀₀ of 0.357, confirming substantial bacterial growth and successful culture revival from cryopreserved stocks.